


# Guide: Setting up the virtual environment to simulate UR robots



Copyright© Universal Robots

All rights reserved. Published 2020

www.universal-robots.com

No part of this publication may be stored in a retrieval system, transmitted, or reproduced in anyway, including but not limited to photocopy, photographic, magnetic or other record, without the prior agreement and written permission of Universal Robots.

This is a registered trademark of Universal Robots. Universal Robots trademarks may not be used in conjunction with any product or service that is not associated with Universal Robots in any manner that is likely to cause confusion amongst customers, or in any manner that disparages or discredits Universal Robots.

Universal Robots shall not be liable for errors contained herein, or for incidental or consequential damages in connection with the use of this material. The contents of this document may be changed without notice.

# Contents

Premise:	1
Part 1: Download of Oracle VM	1
Part 2: Download of URsim Virtual Disk Image	2
Part 3: Set up of the Virtual environment	4
Part 4: Settings for correct execution of the VDI	8
Part 5: VDI first launch	11
Part 6: Launch Polyscope	15
Part 7: Change Keyboard Layout	16
Part 8: Set Up of Network Adapter	17

### **Premise:**

URSim is a simulation software that is used for offline programming and simulation of robot programs and is made for the Linux operating system. This guide shows how to install, setup and run URSim in the MS WINDOWS environment. Please note that the guide is based on VM VirtualBox 6.1.6 and Polyscope 5.8; if you are using previous versions a few different settings will be required.

Universal Robots makes available for users a virtual disk image of a Linux operating system containing the software that simulates Polyscope and the physical robots, called URSim.

To run the simulator in another operating system, a virtual machine is needed.

# Part 1: Download of Oracle VM.

First, in order to be able to launch a VDI (virtual desktop interface) we need software that is capable to do so. Many options are available on the market; for this tutorial Oracle VM VirtualBox is chosen as it is currently free of charge. (<u>https://www.oracle.com/technetwork/server-storage/virtualbox/overview/index.html</u>)

#### Step 1

Go to (<u>https://www.oracle.com/virtualization/technologies/vm/downloads/virtualbox-downloads.html</u>) and download the latest Oracle VM VirtualBox installer.

#### Oracle VM VirtualBox Base Packages - 6.1.6

Freely available for Windows, Mac OS X, Linux and Solaris x86 platforms under GPLv2:

Platform	64-bit
Windows	Uindows Installer
Mac OS X	dmg Image
Solaris 10 5/08 and later or Solaris 11	Solaris Package

Linux Platforms

#### Step 2

Once the download has finished, launch the downloaded installer and follow the suggested procedure.

# Part 2: Download of URsim Virtual Disk Image


Once Oracle VM VirtualBox has been installed, we need to download the Virtual Disk Image that contains the software provided by Universal Robots to simulate the robots.

To do so, go to <u>https://www.universal-robots.com/download</u>, and select the following options from the menu:

- eSeries (or CB3 if you need the previous robot version)
- Software
- Offline Simulator
- Non-Linux
- URSim for non-Linux version 5.8 (or newer)

Select robot type		e-Seri	es 🗸
Label on the robot arm and control box will sh	how what type it is.		
Select type of download fo	or e-Series		
Search for article "How to use this Support si		( <mark>Softw</mark>	are V
Select type of software			
Robot software should ONLY be installed by t The software is not backwards compatible.	trained personnel.	Offline Sir	nulator v
Software is covered under End User Software under software selection.	e Licence Agreement. Please find thi	s	
Select operating system fo	· ·		
URSim is a simulation software intended for o both robot programs and manual movement		of	
NOTE: Not all functions works compared to a - Emergency stop can not be used - Input IO state can not be set - Paths are perfect - Collisions with self or with surrounding obje - Force mode will not work		Noni	Linux v

Then, at the bottom of the page, click on URSim\_VIRTUAL-5.8.0.10253.rar to start your download.




Once the download has finished, extract the downloaded .rar file to the directory of choice.

· · · ·	iis PC > OSDisk (C:) > Users > pica >	Downloads						ٽ ×	,○ Search
🖡 Downloads  🖈 🐴	Name		Date modified	Type		Size			
🖥 Documents 🚿	V Today (1)								
E Pictures 🛛 🖈	@ URSim_VIRTUAL-5.8.0.10253		17/04/2020.00.42	0.40.5%		822.227 KB			
eSeries	✓ Last month (14)	Open							
ML CB3	python-3.8.2		arm Community Edition		on	25.861 KB			
ML eseries	ISS project	<ul> <li>Move to Onel</li> </ul>	Drive	1	on	23.001 ND			
UR_DOC	Diagnostics_Logs-OLK-UTC.2020	7-Zip		>	Apri				
OneDrive - Univer-	e_Series_Brochure_IT_WEB	CRC SHA		>	Apri			>	
Desktop	📴 7z1900 (1)	Edit with Note	pad++		Estrai i files	i			
Documents	🌆 7z1900	Open with		5	Estrai qui	RSim VIRTUAL-	0.0.4025218		
Pictures	Gestione Piani di manutenzione				Verifica l'ar	-	5.8.0.10253\		
	KMLRPC	Give access to		>		ll'archivio			
This PC	python-2.7.amd64	Restore previo	ous versions			ed invia per ema	il		
3D Objects	URCap-getControllerTime-master Robotics in UR Robots March 20	Send to		>			AL-5.8.0.10253.rar.7z"		
Desktop	JavaSetup8u241 (1)	Cut			00 0	-	JAL-5.8.0.10253.rar.7z" ed invia per email		
Documents	JavaSetup8u241 (1)	Сору			Aggiungi a	"URSim_VIRTUA	AL-5.8.0.10253.rar.zip"		
Downloads	Gestione Piani di manutenzione	Create shortcu	ıt		Comprimi	in "URSim_VIRTL	JAL-5.8.0.10253.rar.zip" ed invia per email		
Music	✓ Earlier this year (7)	Delete		Ĭ					
Pictures	URLAWI	Rename				18 KB			
🖞 Videos		Properties			sed (zipped)	6.683 KB			
) OSDisk (C:)	🖡 dist		21/02/2020 14:48		sed (zipped)	30.152 KB			
Common (R:)	Jubuntu-16.04.6-desktop-i386		21/02/2020 09:12	Disc Imag	ge File	1.638.400 KB			
🙀 Faelles (X:)	Starter+Package+SDK+1.9.0		19/02/2020 15:11	Compres	sed (zipped)	10.101.081 KB			
Network Y	ELECTROLUX		24/02/2020 14:23	File folde	ar.				

(If needed, you can find a RAR extractor at https://www.7-zip.org/download.html)

# If the procedure has been executed correctly, the directory to which you extracted the files should look like this:



# Part 3: Set up of the Virtual environment

Once all the needed software and files have been downloaded, we must set up Oracle VM VirtualBox to launch the URSim VDI.

#### **Step 1:** Launch Oracle VM VirtualBox.

🜍 Oracle VM VirtualBox Manager File Machine Help	_	□ ×
Tools	E Contraction of the second se	
Gatto O Powered Off	Welcome to VirtualBox! The left part of application window contains global tools and lists all virtual machines and virtual machine groups on your computer, You can inport, add and create new VMs using corresponding toolbar buttons. You can pool a tools of currently	
Corso_ros_pica	computer. You can import, add and create new VMs using corresponding toolbar buttons. You can popup a tools of currently selected element using corresponding element button. You can press the F1 key to get instant help, or visit www.vitualbox.org for more information and latest news.	
PythonCompiler Powered Off		
UrCap Dev env () Powered Off		

**Step 2:** Add a new VBoxImage by clicking on the 'New' button.

👽 Oracle VM VirtualBox Manager	
File Machine Help	
Tools	New Settings Discard Show
Gatto Powered Off	General           Name:         URSIM_5.8           Operating System:         Ubuntu (64-bit)
Corso_ros_pica	System Base Memory: 1024 MB Boot Order: Floppy, Optical, Hard Disk
PythonCompiler O Powered Off	Acceleration: VT-x/AMD-V, Nested Paging, KVM Paravirtualization Display
UrCap Dev env	Video Memory:     16 MB       Graphics Controller:     VBoxVGA       Remote Desktop Server:     Disabled       Recording:     Disabled
₩ URSIM_5.8 Running	Storage Controller: IDE IDE Secondary Master: [Optical Drive] Empty Controller: SATA SATA Port 0: URSim_VIRTUAL-5.8.0.10253.vmdk (Normal, 10,00 GB)
	Audio           Host Driver:         Windows DirectSound           Controller:         ICH AC97

#### Step 3: Select Name and Operating System

		?	$\times$
Create Virtual I	Machine		
Name and	operating system		
machine and sel	ect the type of operating system you inte	end to install o	on it.
Name:	URSIM		
Machine Folder:	C:\Users\pica\VirtualBox VMs		~
Type:	Linux	-	64
Version:	Ubuntu (64-bit)	-	-
	Expert Mode Next	Cano	cel
	Name and of Please choose a machine and sel The name you of machine. Name: Machine Folder: Type:	machine and select the type of operating system you inte The name you choose will be used throughout VirtualBox machine. Name: URSIM Machine Folder: C:\Users\pica\VirtualBox VMs Type: Linux Version: Ubuntu (64-bit)	Create Virtual Machine Name and operating system Please choose a descriptive name and destination folder for the new vi machine and select the type of operating system you intend to install of The name you choose will be used throughout VirtualBox to identify th machine. Name: URSIM C:\Users\pica\VirtualBox VMs Type: Linux Version: Ubuntu (64-bit)

While the name can be chosen freely, in this case URSIM, the Type entry (Linux) and the Version (Ubuntu 64 Bit) must be chosen as shown.

V

Once finished, click 'Next'. Then select a sufficient amount of RAM memory:


2

		1	
←	Create Virtual Machine		
	Memory size		
	Select the amount of memory (RAM) in megabytes to be allo virtual machine.	cated to the	
	The recommended memory size is <b>1024</b> MB.		_
		1024 🗘	MB
	4 MB 8192 MB		_
	Next	Cance	
	NEAL	Curice	

At least 768 MB should be selected; if possible, reserve more than that, specifying a value equal to a power of 2 (1024, 2048 etc.).

Click 'Next'.


In the next window, select 'Use an existing virtual hard disk file'.




Click on the folder with the green arrow on the right, and in the new window select 'Add':

<u> aa</u> - Hard Disk Selector		?	$\times$				
Medium							
Add Refresh							
Name	Virtual Size	Actual Size	^				
✓ Attached							
corso_ros_pica.vdi	15,00 GB	12,01 GB					
Gatto.vdi	30,00 GB	11,43 GB					
PythonCompiler.vdi	10,00 GB	9,14 GB					
URCaps Starter Package.vmdk	20.00 GB	19.04 GB	~				

Navigate to the folder where you extracted the URSIM file downloaded from the UR site and select the first element:



Click 'Open', then 'Choose' and then 'Create':




If everything has been done correctly in the Oracle main page on the left side, you should have a new entry with the selected name:

🜍 Oracle VM VirtualBox Manager		- 🗆 ×
File Machine Help		
Tools	New Settings Discard Start	
Gatto	General Name: URSIM Ubuntu (64-bit) Settings File Location: C:\Uesespica\VirtualBox VMs\URSIM	Preview
Powered Off	System Base Memory: 1024 MB	URSIM
PythonCompiler Powered Off	Boot Order: Floppy, Optical, Hard Disk Acceleration: VT-x/AMD-V, Nested Paging, KVM Paravirtualization	
UrCap Dev env Opwered Off	Display	
VIRSIM O Powered Off	Video Memory: 16 MB Graphics Controller: VMSVGA Remote Desktop Server: Disabled Recording: Disabled	
~	Storage	
	Controller: IDE IDE Secondary Master: [Optical Drive] Empty Controller: SATA SATA Port 0: URSim_VIRTUAL-5.8.0.10253.vmdk (Normal, 10,00 GB)	
	🖗 Audio	
	Host Driver: Windows DirectSound Controller: ICH AC97	
	📑 Network	
	Adapter 1: Intel PRO/1000 MT Desktop (NAT)	
	🖉 USB	
	USB Controller: OHCI Device Filters: 0 (0 active)	
	Shared folders	

# Part 4: Settings for correct execution of the VDI

In order to correctly run the VDI we must first set some preferences; to do so, select the VDI installed before and click on 'Settings':



The following window is opened:

😟 ur	SIM - Settings							?	×
	General	Gener	al						
	System	Basic	Advanced	Description	Disk Encryption				
	Display	Name:	URSIM						
$\mathbf{\Sigma}$	Storage	Type:	Linux					•	
	Audio	Version:	Ubuntu (64-t	pit)				•	
5	Network								
	Serial Ports								
Ď	USB								
	Shared Folders								
	User Interface								
							ОК	Ca	incel

Click on 'Display' on the left side of the window to access the following window:

🙆 UR	SIM_5.8 - Settings		?	×
	General	Display		
	System	Screen Remote Display Recording		
	Display	Video Memory:	16 MB	-
$\mathbf{D}$	Storage	0 MB 128 MB 128 MB	1	* *
	Audio	1	8	
•	Network	Scale Factor: All Monitors - Min Max	100%	•
	Serial Ports	Graphics Controller: VBoxVGA •		
Ď	USB	Acceleration: Enable 3D Acceleration		
	Shared Folders			
	User Interface			
		Invalid settings detected OK	Car	ıcel

#### Change the Graphics Controller to VBoxSVGA:


🙆 UR	SIM_5.8 - Settings	;	?	×
	General	Display		
	System	Screen Remote Display Recording		
	Display	Video Memory:	16 MB	* *
	Storage	0 MB 128 MB		\$
	Audio	1	8	
-	Network	Scale Factor: All Monitors 🔹 🖣	x 100%	÷
	Serial Ports	Graphics Controller: VBoxSVGA -		
	USB	Acceleration: Enable 3D Acceleration		
	Shared Folders			
	User Interface			
		Invalid settings detected 🔥 OK	Can	ncel

Please ignore the warning "Invalid settings detected". Some users experienced issues with the screen size if the default "VMSVGA" option was choosen.

Optionally, we need to set up a shared folder between the Virtual System and our Windows System.

This part is necessary only if you plan to transfer files from the virtual machine to Windows environment or vice versa, and can be done at a later time (Jump to chapter 5, if not needed).

On the left side of the window, select 'Shared Folder' and click on the folder with a green cross on the right side of the window:



Specify the options as follows:

- **Folder path**: Specify the path to a folder in your Windows system: This is the Windows folder where you will put files you need to transfer to the simulator in the Linux system.
- Folder Name: Choose a desired folder name
- Check the Auto Mount
- Mount Point: Specify the mount point in the Linux system from where you will be able to access the content of the shared folder, remember to start with a "~" (the "~" can be typed as "Alt"+"126" on the keypad)


Example:

😟 Edit Share		?	$\times$		
Folder Path:	C:\Users\pica\VirtualBox VMs\New folde	h	~		
Folder Name:	lder Name: New_folder				
	Read-only Auto-mount				
Mount point:	~/sharedvm				
	OK	Cano	cel		


Click 'OK' and again click 'OK' to close the Settings window, we are now ready to launch our VDI.

# Part 5: VDI first launch

In Oracle, select your new VDI and click on 'Start':



You should now be inside your new Linux OS:



If the screen is small size, expand the window clicking on the right top square:




And drag (one click and hold) the URSim UR16 and Programs UR16 icons to correctly visualize all the icons on the desktop:



**The following part is necessary only if you defined a shared folder in the previous chapter.** (Else jump to chapter 6).

In order to be able to access the shared folder created in chapter 4, open a terminal and type the following line:

sudo adduser ur vboxsf (see chapter 7 if you need to change the keyboard layout)




Press Enter, close and restart the virtual machine.

Once restarted, the shared folder is located in the directory " /media "


To access it, click on the icon in the bottom left corner, then in the window move your cursor to accessories and click on File Manager PCManFM

	ograms URS
URSIM UR10	Programs UR16
URSim UR16	Programs File Manager PCManFM Con poo Programs UR10
il <sup>®</sup> Accessories	Programs UR16
<ul> <li>System Tools</li> <li>Preferences</li> <li>Run</li> </ul>	, III Programs UR3 , III Programs UR5 III URSIM UR10 III URSIM UR16
Logout	URSim UR3


The following window will be opened; click twice on the up-arrow shown in figure...



...to reach the following directory, then double click on media:



#### Here you will find your shared directory:




Double click to see its contents. Your machine should now be ready for use, the screen should be of the correct size, and you should be able to access the shared folder between the Virtual Machine and your Windows machine to easily transfer files between the two.

## Part 6: Launch Polyscope

Launching Polyscope and the whole robot simulation is very easy.


You can choose to launch a simulation of one of the four different types of UR Cobots: UR3e, UR5e, UR10e and UR16e.

Depending on your choice from the desktop of your Virtual Machine, double click on the desired icon:



The folders 'Programs URx' contain .urp programs, .scripts and .urcap used by the simulator. There you will find all the software created in the simulator and this is where you will put all the files (from a real robot or received by e-mail) that you need to run in the robot simulator.


Once you launch one of the available simulations, a window will appear as if you are holding the Teach Pendant of the real robot:



# Part 7: Change Keyboard Layout

This part is necessary only if you need to change the keyboard layout (else jump to chapter 8).

By default the keyboard layout of the virtual machine is the US one, if you have a different layout you can change it by clicking on the US symbol on the bottom right corner of the screen:



In the following window, click on "add" to add a new keyboard layout:

	Keyboard Layo	ut Handler			- + ×		
Keyboard Model		Advanced set	xkbmap	Options			
					č		
pc105		Do not res					
		🗌 Keep syste	em layou	uts			
Keyboard Layouts	-				Add Keybo	oard Layout	- + >
Flag Layout Variant	+ Add	Per Window S	Flag 🛩	Layout	Description		*
us	- Remove		• 6		Afghani		
		Remember		al	Albanian		
	∧ Up	Show Layout	<b>T</b> io <b>T</b>	et	Amharic		
	Y Down		• @	ага	Arabic		
				ma	Arabic (Morocco)		
		🖾 🔾 Cust		sy	Arabic (Syria)		
		US 💿 Text	۰ <b></b>	am	Armenian		
Change Layout Option			<u>همه</u>	az	Azerbaijani		
		Panel Icon Siz	۰ <b></b> 📕	ml	Bambara		
shift_caps_toggle		01 02	۰ 💻	by	Belarusian		
		01 02	۰ <b>🔳</b>	be	Belgian		
			› 💽	bd	Bengali		
			۰ 📐	ba	Bosnian		
			۰ 🚃	brai	Braille		
			۰ 💻	bg	Bulgarian		
				mm	Burmese		
			2	ad	Catalan		
			<ul> <li>•</li> </ul>	cn	Chinese		
			• ==	hr	Croatian		
							X Cancel ✓ OK

Clicking on the Up / Down buttons you Will change their priority, the one on top will be the one used as default.

# Part 8: Set Up of Network Adapter


ADVANCED USERS - This part is necessary only if you plan to establish a network connection with the WINDOWS environment or to external devices (else, stop here)

In settings \ network, depending on the hardware used and the need to comuncate with a software in the host or an external device, you should chose to connect to "Host Only Adapter" or "adapter with bridge".

Your WINDOWS environment and your URSIM virtual machine Will now have two differents IP addresses that can be used to comunicate.

The IP address of your WINDOWS environment used to comunicate FROM the Virtual Machine TO the Windows environment can be retrieved as follows:

- In the bottom left, in the search text box type "cmd" and the start the command prompt:




- Once the command prompt is opened type the command "ipconfig" and press enter:
- Search for the "Ethernet Adapter Virtual Box Host Only network", the IP address shown in here is the one that you will need to use when comunicating from the virtual machine to the Windows environment:




In the same way, to know the IP address of your Virtual Machine to be used to communicate to it from the WINDOWS environment follow this procedure:


- Start your virtual machine
- In the bottom right corner click on the icon highlighted in the figure, then select System Tools and the UXTerm, as shown before;
- The Linux terminal Will open, now type the command "ifconfig":



- The IP address of your virtual machine is the one shown inside eth0:



To verify the connection is properly configured from the Windows Command prompt type the following command "ping #IPADDRESSOFVIRTUALMACHINE" and verify you have a response:



Viceversa from the Virtual Machine terminal type the command "ping <IPADDRESSOFWINDOWS>":

🛣 ur@ursim: ~	- + ×
ur@ursim:"\$ ping 192,168,56,102 PING 192,168,56,102 (192,168,56,102) 56(84) bytes of data. 64 bytes from 192,168,56,102: icmp_seq=1 ttl=128 time=0,333 ms 64 bytes from 192,168,56,102: icmp_seq=2 ttl=128 time=0,236 ms 64 bytes from 192,168,56,102: icmp_seq=4 ttl=128 time=0,387 ms 64 bytes from 192,168,56,102: icmp_seq=5 ttl=128 time=0,393 ms 64 bytes from 192,168,56,102: icmp_seq=5 ttl=128 time=0,393 ms 64 bytes from 192,168,56,102: icmp_seq=5 ttl=128 time=0,393 ms 64 bytes from 192,168,56,102: icmp_seq=6 ttl=128 time=0,504 ms ^C 192,168,56,102 ping statistics 6 packets transmitted. 6 received. 0% packet loss, time 4997ms rtt min/avg/max/mdev = 0,236/0,364/0,504/0,083 ms ur@ursim:"\$	